

Alert Level: **INFORMATIONAL**

Instruments: All Profilers and Samplers

Software Version: N/A

Hardware Version: N/A

**Subject:** Recommended USB to RS-232 Adaptor for Connecting a Computer to McLane Instruments

**Summary:** A USB to RS-232 adaptor cable is required to connect a McLane Profiler or Sampler to a computer without a built-in serial port.

**Technical Details:** Computers with no built-in serial port need a USB to RS-232 adaptor cable to connect to McLane instruments. McLane recommends the DS\_US232R-10(R-100/500) adaptor cable manufactured by FTDI Ltd. This adaptor cable is available from McLane or other commercial suppliers.

**Actions:** Contact McLane for a USB to RS-232 adaptor cable or order the adaptor cable from a commercial supplier such as Digikey (p/n 768-1013-ND).



**FTDI Ltd**

**DS\_US232R-10(R-100/500)**

**USB to RS232 Adapter cable**