

**Alert Level: INFORMATIONAL**

**Instruments: WTS-LV (all models)**

**Software Versions: N/A**

**Hardware Version: N/A**


**Subject: MCBH Bulkhead Connector Upgrade**

**Summary:** Replacing any 'legacy' (Teledyne Impulse AG and XSG series) cables and connectors will require an upgrade to Subconn MCBH cables and connectors.

**Technical Details:** Legacy' Impulse AG and XSG cable and connector. Customers who need to replace an old-style cable or connector must upgrade to Subconn MCBH cables and connectors. The full MCBH upgrade for WTS-LV connectors includes the following:

- New pump cable
- New COM cable
- New COM dummy plug
- New MCBH connector with Ti adapter for pump to end cap connection
- New Trigger plug (if Trigger option is installed)

Figure 1 shows the parts for the upgrade. Contact [mclane@mclanelabs.com](mailto:mclane@mclanelabs.com) for upgrade pricing.



*Figure 1: MCBH cable and connector upgrade*