WTS-LV Filter Types and Recommended Pump Heads

The tables shown next provide a list of available pump heads and recommendations based on compatible filter type and flow rate.

Pump Size	LV04	*LV08	LV30
Flow Rate Range (L/min)	1–5	4–10	12–30

* LVO8 is standard on the WTS-LV

Recommended Pump Heads

Filter Type	Pore Size in µm	Max Flow Rate (L/min)	LV04	LV08	LV30
Polycarbonate and Polyester Membrane (Nuclepore®, Millipore®)	0.2	2	•		
Polycarbonate and Polyester Membrane (Nuclepore®, Millipore®)	0.4	3	•		
Polycarbonate and Polyester Membrane (Nuclepore®, Millipore®)	0.8	6	•	•	
Polycarbonate and Polyester Membrane (Nuclepore®, Millipore®)	1.0	8	•	•	
Glass Fiber (GF/F®, QMA®)	0.8	5 (LV04 pump) 7 (LV08 pump)	•	•	
	1.7	5 (LV04 pump) 10 (LV08 pump)	•	•	
	5.0	5 (LV04 pump) 10 (LV08 pump) 20 (LV30 pump)	•	•	•
Mesh (Nytex®)	60.0 (+)	5 (LV04 pump) 10 (LV08 pump) 30 (LV30 pump)	•	•	•